



Department/Speciality: \_\_\_\_\_

Ward/Bed: \_\_\_\_\_ Clinic: \_\_\_\_\_ Class: \_\_\_\_\_

Biochemical Genetics
Patient's name label (For Downtime Use) Name: MRN: Account Number: Date of Birth: Sex: M / F (Circle one)

Patient Type <input type="checkbox"/> Gynae <input type="checkbox"/> Obst <input type="checkbox"/> Neo <input type="checkbox"/> Medicine <input type="checkbox"/> Surgery	Relevant History/ Findings/ Treatment	Laboratory Barcode
Clinical Diagnosis		<b>For Laboratory Use Only</b>
Name & Signature of Requesting Doctor Pager / Contact no. ( <i>indicate if urgent</i> ) Name of Consultant I/C Date	Type of Specimen - Please use one form per specimen type <input type="checkbox"/> Urine, Random <input type="checkbox"/> Others <input type="checkbox"/> Blood (LiHep) Specify _____ <input type="checkbox"/> CSF Date & Time Specimen Taken Date _____ Time _____ am / pm	

Please tick ( ) appropriate boxes below

PLEASE NOTE THAT TO REQUEST URGENT ANALYSIS YOU MUST CALL THE LAB AT EXT 8728/5049.

WRITING ON THE FORM IS NOT SUFFICIENT

URINE TESTS (Fresh random sample, no preservative.)		
BC6010U	<input type="checkbox"/>	Amino Acids, urine 3ml fresh urine, no preservative. Please send with a separate ice pack.
BC6090	<input type="checkbox"/>	GAG electrophoresis, urine 5 to 10ml fresh urine (3 ml is the absolute minimum), no preservative. Please send with a separate ice pack. For diagnosis of Mucopolysaccharide disorders.
BC6180U	<input type="checkbox"/>	Methylmalonic Acid Quantitative, urine 3 to 10ml fresh urine, no preservative. Please send with a separate ice pack.
BC6150U	<input type="checkbox"/>	Organic Acids, urine 3 to 10ml fresh urine, no preservative. Please send with a separate ice pack.
BC6160U	<input type="checkbox"/>	Orotic Acid Quantitative, urine 3 to 10ml fresh urine, no preservative. Please send with a separate ice pack.

BLOOD TESTS (Lithium heparin tube)		
BC6010	<input type="checkbox"/>	Amino Acids - Full Quantitation 1ml Lithium Heparin Tube. Send straight to lab as ideally should be separated within 1hr of sampling. Please send with a separate ice pack.
BC6020	<input type="checkbox"/>	Amino Acids - Part Analysis 1ml Lithium Heparin Tube (send straight to lab as ideally should be separated within 1 hr of sampling). Please send with a separate ice pack.
BC6050	<input type="checkbox"/>	Gal-1-Phosphate-Uridyl Transferase 0.5ml Lithium Heparin Tube. Please send with a separate ice pack.
BC6461	<input type="checkbox"/>	Sterol Analysis 1ml Lithium Heparin Tube, keep away from direct sunlight. Send straight to lab. Please send with a separate ice pack. Includes 7 dehydrocholesterol, cholestanol, phytosterols.
BC6462	<input type="checkbox"/>	VLCFA, Phytanic Acid & Pristanic Acid, plasma / Peroxisomal Disorders Panel 1ml Lithium Heparin Tube. Send straight to lab. Please send with a separate ice pack.
BC6070	<input type="checkbox"/>	Acylcarnitines, plasma 1ml Lithium Heparin Tube. Please send in a separate ice pack.

CSF TESTS (plain sterile container, no preservative)		
BC6010F	<input type="checkbox"/>	Amino Acids, CSF 200µl (no preservative) Pls send with a paired plasma sample (collected within an hour) as interpretation, particularly in the case of nonketotic hyperglycinaemia (NKHG) requires comparison with plasma concentrations. Please send with a separate ice pack.

 Others: pls specify \_\_\_\_\_